



## FIRST DRAW & FLUSH LEAD TEST SAMPLING INSTRUCTIONS

### PLEASE READ ALL INSTRUCTIONS COMPLETELY BEFORE COLLECTING YOUR WATER SAMPLE

1. First Draw sample: Choose a collection time when the water has not been used for at least six (6) hours. Before you use any water from the tap, fill the sample (250 mL) vial labeled "First Draw" as soon as you turn on the tap.
2. Flush sample: Run the water from the tap for 3-5 minutes and then fill the sample (250 mL) vial labeled "Flush".
3. Complete the lab form, **making sure to indicate the collection date and time you sampled**, type of water sampled (treated/untreated) and sample identification which identifies where you collected your sample (i.e., kitchen sink, well head, etc.). Please be sure to retain a copy of the completed lab form prior to shipment. Return the completed paperwork and water samples to the shipping container.
4. It is not necessary to ship by an overnight delivery service. Ground, or Saver delivery is enough. Cost of shipping samples to the laboratory is at your expense.

Testing Location: **AIH Laboratory**  
**2556 W Woodland Dr**  
**Anaheim, CA 92801**

5. PLEASE NOTE: Testing is performed during business days based on the turnaround selected from receipt of samples. Results will be emailed upon completion.
6. You may use a clean plastic/glass container to collect.



**NOTE: The lab reserves the right to reject water samples that appear to be excessively contaminated as they may cause potential damage to laboratory equipment.**

If you have any questions regarding this First Draw & Flush Lead Test or any other services we provide, please contact a Technical Service Representative at AIH Laboratory at 1-562-860-2201 for assistance.

Thank you for choosing AIH Laboratory for your water testing services!

**AIH Laboratory**  
Address: 2556 W Woodland Dr  
Anaheim, CA 92801  
Phone: 1-562-860-2201  
Email: [frontdesk@aihlab.com](mailto:frontdesk@aihlab.com)  
Web: [www.aihlab.com](http://www.aihlab.com)